ISP|80 Unit

with 5 equipment rails



I-PA2227

PRODUCT DESCRIPTION

| I KODOOI DEGOK | 11014 |
|-------------------------------|---------------------------------------|
| Material | Stainless steel |
| Bottle Holder | Stainless steel, bottle bar |
| Number of hooks | 6 |
| Arrangement of the hooks | In line with each other |
| Type of height adjustment | Screw height adjustment |
| Maximum length/height | 2,200 |
| Minimum length/height | 1,690 |
| Inner tube diameter | Ø 25 mm |
| Outer tube diameter | Ø 38 mm |
| Length of outer tube | 1,420 |
| Base leg profile | 45 x 25 mm |
| Base weight | Weighted |
| Base diameter | Ø 680 mm |
| Castors | Twin castors |
| Castor diameter | Ø 75 mm |
| Number of castors | 5 with brakes and elect. conductivity |
| Total load capacity (kg) | 25 |
| Load capacity per system (kg) | 12 |
| Load capacity per hook (kg) | 2 |
| | |

PACKAGING UNIT: UNIT

| Width Box1 | 0,31 m |
|-------------|---------|
| Height Box1 | 0,21 m |
| Length Box1 | 1,59 m |
| Weight Box1 | 25,7 kg |



Best hygienic properties stainless steel (1.4301)







I-PA2227

ISP|80 Unit









Area of application

comfortable height adjustment through a screwheight adjustment with an integrated, surface-gently mold clamping element

The outer tube offers enough space to attach various devices at a userfriendly operating height



Technology

All plastic parts are lightfast

Stainless steel tubes according to DIN EN 10217-7

Material ultrasonic and eddy current tested (according to EN 10246-9/3)

The integrated pull-out limiter prevents the adjustable tubes from separating

Safety bottle cross according to DIN ISO 15375

great running characteristic through high value 80 mm Mono castors

a sealing cap prevents the IV-pole from hygienical senstive liquids

Reinforced bottle holder allows hanging 5 liters of rinsing solutions



Safety

The base legs are sealed with plastic plugs (no risk of injury)

Low risk of pinching when retracting and extending the IV-pole thanks to the built-in safety distance

an integrated base-weights facilitates a safer attachment of various pumps and other medical equipment

No risk of injury thanks to the rounded shape of the bottle hooks

The geometry of the hooks prevents the infusion bottle from falling out

