

Easy online
configuration:

[https://configurator.
provita.de/endoprotect](https://configurator.provita.de/endoprotect)



provita[®]
medical

First thermo-
disinfectable
protector for validated
reprocessing up to 95 °C!
(for 500 and 700 mm
versions)



ENDO PROTECT

PROTECTED & SAFE

Transport and storage
for endoscopes and fibroscopes


ENDO PROTECT

Transport and storage of:

- video endoscopes
- endoscopes with supply line
- fibroscopes



Technical specifications
at provita.de



• SUSPENDED STORAGE as per RKI REQUIREMENTS**

• IMMEDIATE AVAILABILITY

of the reprocessed endoscope
(e.g. for emergency intubation)

• CONFIGURABLE

suitable for all plug brands, lengths and supply lines

• SEPARATE DUST PROTECTOR

available for the handle

- Available as a **COMPLETE SYSTEM**,
possible protector installation options:

- » infusion stand*
- » docking cart
- » wall rail
- » cabinet installation
- » operating table
- » ceiling supply unit
- » multi-purpose cart

• MDR-COMPLIANT

* For mobile use, it must be ensured that the unit does not tip over and that it is mounted in a stable position, etc.

** Robert Koch Institute



immediately
ready for use

suitable for
all brands

thermo-
disinfectable and
wipe-disinfectable
protector




mechanical
protection of
the handle
during storage
and transport



safe storage and
transport of
the supply plug

protection of
the distal end



THE HYGIENE CONCEPT

The protector enables flexible and safe in-hospital transport of the endoscope from storage to examination. The protector ensures immediate availability of the reprocessed endoscope on the emergency trolley or IV pole. This allows for a quick response in emergencies.





Figure shows thermo-disinfectable protector (versions 500 and 700 mm)

Mobile systems

MOBILE BASE UNITS



for 1 to 2 protectors

T-W02000



for 3 to 4 protectors

T-W04000



PROTECTOR
500 mm

E-C11109



PROTECTOR
700 mm

E-C11119



PROTECTOR
1,200 mm

E-W11129

DISPOSABLE DUST PROTECTION



reprocessed (green)


Z-EC0010



contaminated (red)

Z-EC0020

EQUIPMENT CARRIER SUPPORT



for protector in the WD

Z-EC0002

PUSH HANDLE



Z-Z00002

HANDLE SUPPORTS

supplied pre-assembled in the protector



Z-EH-001



Z-EH-002



Z-EH-003



Z-EH-004



Z-EH-005



MANY other options and
dimensions at provita.de

 HANDLE MOUNTS



small

included with the protector Z-HA0001



medium

included with the protector Z-HA0005



extra-high version

Z-HA0010

 PLUG HOLDER
pre-assembled on the protector



conical

Z-ES-000



cylindrical

Z-ES-011



rectangular

Z-ES-014



VARIOUS other plug holders
and dimensions at provita.de

 SWING GUARD



10 mm

Z-EP-001




23 mm

Z-EP-002




high version

Z-EP-003



Stationary systems

STATIONARY SIMPLE-FIX SYSTEMS



WITH MAGNET



for handle mounts and plug holders,
2 kg load capacity

E-B-S001



for swing guard,
0.5 kg load capacity

E-B-S002



for handle mounts and plug holders,
2 kg load capacity

E-B-S003

WITH SCREWS

WITH UNIVERSAL CLAMP



vertical installation

E-B-S005



horizontal installation

E-B-S004

HANDLE SUPPORTS (8 MM) for use in handle mounts (see below)



Z-EHS001



Z-EHS002



Z-EHS003



Z-EHS004



Z-EHS005

HANDLE MOUNTS for optional handle supports (see above)



left-hand version

Z-EAS010



right-hand version

Z-EAS011



MANY other options and
dimensions at provita.de

SWING GUARD



10 mm

Z-EPS001



23 mm

Z-EPS002



high version

Z-EPS003

HANDLE MOUNTS



conical

20/16 mm Z-EAS001


24/20 mm Z-EAS002

32/28 mm Z-EAS003



universal

Z-EAS004



ENT 12 mm

Z-EAS005



Universal, with mount for supply cable, left-hand version

Z-EAS006



Universal, with mount for supply cable, right-hand version


Z-EAS007

PLUG HOLDER



conical

Z-ESS000



cylindrical

Z-ESS011



rectangular

Z-ESS014



VARIOUS other plug holders
and dimensions at provita.de

PULLEY



Z-Z00360




More options at provita.de

Handle mounts and plug holders are available for many endoscopes:



- + Philips
- + GE
- + Siemens
- + Canon
- + Fujifilm Sonosite
- + Aloka Hitachi
- + Karl Storz




Please contact us –
we can develop a custom
handle mount and plug
holder for you.



Easy online configuration:
<https://configurator.provita.de/endoprotect>

Any further questions?
We are happy to help you personally.

Just call us at: +49 (0) 21 93 / 51 05 - 0
Or write to us: info@provita.de





Have you seen our other catalogues?



Version: 2023

provita medical gmbh & co. kg
Auf der Huhfuhr 8
42929 Wermelskirchen
Germany

Phone +49 (0) 21 93 / 51 05 - 0
Fax +49 (0) 21 93 / 51 05 - 269

info@provita.de
www.provita.de

