



## Patient transfer

**SAM TRANSFER CHAIR** 

Your • in service and hygiene

#### SAM



## **Mobility in** health and care

The modern transfer chair for many applications for indoor and outdoor use



S-SAM100 Indooruse

Central lockable foot brake,

steel and powder coated ......p. 4





S-SAM001 Indoor and outdoor use

Dead man brake, washable, 





S-SAM004 Indoor and outdoor use

Similar to 001, but with an advertising set ...... p. 7





S-SAM003 Indoor use

Outpatient version with leg support ...... p. 8

#### This applies to all models





is possible thanks to folding armrests

### can be loaded up to 150 kg



The safe automatic footrest (glass fibre reinforced) enables



Make it your own chair with



seat and back parts made of polyurethane are



#### **S-SAM100**

## The convenient indoor aid

with central lockable foot brake made of powder coated steel





#### SAM Indoor

**Frame material** steel tubing  $\emptyset$  32 × 2 mm powder coated

Material for armrests, backrests and seats PU foam


**Dimensions** approx. 960 × 710 × 990 mm (L × W × H)

Weight 21.0 kg

S-SAM100








Stacking dimension =  $[(n-1) \times 305 \text{ mm}] + 970 \text{ mm}$ n = number of chairs



Space saving and stackable



#### S-SAM001/S-SAM004

## The washable, low maintenance chair

Dead man brake, washable, maintenance-free PU tyres





### SAM with dead man brake

**Frame material** Stainless steel tubing 18 – 10, Ø  $32 \times 2$  mm electropolished

Material for armrests, backrests and seats PU foam

**Dimensions** approx. 960 × 710 × 990 mm (L × W × H)

Weight 21.0 kg

S-SAM001



### SAM with advertising set

Frame material Stainless steel tubing 18 – 10, Ø 32  $\times$  2 mm electropolished

Material for armrests, backrests and seats PU foam

**Dimensions** approx. 960 × 710 × 990 mm (L × W × H)

Weight 21.0 kg

**Advertising set comprising** 2 stickers on the rear wheel fork and 2 stickers on the side advertising plate

S-SAM004



Washable



Safe and comfortable

#### **S-SAM003**

# Support in an emergency

Transport chair for outpatient use with variable leg support







## SAM Outpatient version

**Frame material** Stainless steel tubing 18 - 10, Ø  $32 \times 2$  mm electropolished, coated in RAL 3020 - traffic red

Material for armrests, backrests and seats PU foam

**Dimensions** approx.  $960 \times 710 \times 990 \text{ mm } (L \times W \times H)$ 

Weight 23.0 kg

S-SAM003







#### **Accessories**Safety and storage





#### LOCK\*

For anti-theft protection

Z-SAM019

#### **SAFETY BELT**

Fast and safe

Z-SAM002



#### **DEPOSIT SYSTEM\***

Coin

Z-SAM003



#### **DEPOSIT SYSTEM\***

Card



#### **TILT SUPPORT**

Strongly recommended when using the oxygen holder or equipment rail

Z-SAM014



#### **REPLACEMENT KEYCARD**

Stainless steel

Z-SAM010



#### **PARKING STATION**

**Free standing** 

Z-SAM004



#### **PARKING STATION**

For wall mounting





**LUGGAGE HOLDER** 

**Practical and sturdy** 

Z-SAM011



#### **DOCUMENT HOLDER**

**Easy to mount** 

Material aluminium

Inner dimensions 330 × 85 mm

Z-SAM006



**MULTI-PURPOSE BAG** 

**Especially for the transfer chair** 

Z-SAM008



**STANDARD RAIL**\*

Only in conjunction with the tilt support

#### **Accessories**Medical care/hygiene



**HOLDER**\*

For an oxygen bottle with 10 cm diameter

Z-SAM012



#### **WALKING STICK HOLDER**

Suitable for various walking aids

Z-SAM013



#### **INFUSION STAND**

**Especially for the transfer chair** 

Z-SAM001



#### **STAINLESS STEEL PUSH HANDLE**

For S-SAM100



<sup>\*</sup> We strongly recommend using the tilt supports (Z-SAM014) when using the oxygen holder and standard rail.



#### **ADVERTISING STICKER SET**

"Light"

Z-SAM016



#### **ADVERTISING SET**

**Signs** 

Z-SAM017



#### **ADVERTISING SET**

Comprising "Z-SAM017" and "Z-SAM016"























Our general terms and conditions apply. Errors, misprints and changes reserved. © Copyright: provita medical gmbh & co. kg, 2023. No part of this publication may be reprinted or reproduced in any form or by any means without the permission of provita medical.

## Any further questions? Please contact us for personal advice.

Just call us on: +49 21 93 / 51 05 - 0

Or write to us at: info@provita.de

provita medical gmbh & co. kg Auf der Huhfuhr 8 42929 Wermelskirchen Germany

Phone +49 (0) 21 93 / 51 05 - 0 Fax +49 (0) 21 93 / 51 05 - 269

info@provita.de www.provita.de



