

https://configurator. provita.de/tee-care







PROTECTED & SAFE

Transport, storage and reprocessing solution for TEE ultrasound probes





- NO BENDING thanks to suspended storage
- MECHANICAL PROTECTION of the TEE probe: protection against damage and swing guard
- NO CONTAMINATION thanks to dust-protected storage
- CLEAR MARKING probe reprocessed/contaminated
- WIPE DISINFECTION smooth surfaces and good accessibility for wipe disinfection
- IMMEDIATE AVAILABILITY of the reprocessed TEE probe
- CONFIGURABLE suitable for all probe brands



- CABLE PULLEY for optimal & kink-free cable guiding
- Available as a COMPLETE SYSTEM, possible protector installation options:
- » cabinet installation
- » IV pole*
- MDR-COMPLIANT

- » docking cart
- » wall rail

* For mobile use, it must be ensured that the unit does not tip over and that it is mounted in a stable position, etc.






THE MOBILITY CONCEPT

The protector enables flexible and safe in-hospital transport of a TEE probe from storage to examination. The protector ensures immediate availability of the reprocessed probe on an emergency trolley or IV pole. This allows for a quick response in emergencies.

Reprocessed probe – stored in the protector and stowed away in a cabinet



Storage of the reprocessed protector in a cabinet or again directly on an IV pole/ emergency trolley






To the video at provita.de/video/TEE/

Attaching the protector to the IV pole or emergency trolley





Probe suitable for immersion disinfection



Reprocessing of the probe using immersion disinfection and reprocessing of the protector using wipe disinfection*

- 1. Cleaning
- 2. Immersion disinfection
- 3. Leak test
- 4. Repeated immersion disinfection
- 5. Storage in the protector

*required process steps depending on specifications and use






Reprocessing of the probe and protector using wipe disinfection



Attachment options/ variants



Mobile unit wipe disinfection





	T-W05000	T-T05030	T-W05040	T-T05050	T-T05060	T-T05070
Enzyme pre-cleaning					✓	✓
Immersion disinfection		✓		✓	✓	✓
Leak test			✓	✓		✓
Transport	✓	✓	√	√	✓	√



Easy online configuration:

https://configurator.provita.de/tee-care



Stationary for standard equipment rails

Stationary units













	T-W55000	T-T55030	T-W55040	T-T55050	T-T55060	T-T55070
Storage	✓	✓	✓	✓	✓	✓
Immersion disinfection		✓		✓	✓	✓
Leak test			✓	✓		✓
Enzyme pre-cleaning					✓	✓

Stationary units Stationary units

FOR MOUNTING ON STANDARD EQUIPMENT RAILS

Protector for standard equipment rail

 Green insert
 T-W71416


 Red insert
 T-W71426

Plug holder for standard equipment rails

see page 14



10



FOR CABINET MOUNTING

Cabinet unit


for 1 ultrasound probe	T-W06401
for 2 ultrasound probas	T W06402








Easy online configuration: https://configurator.provita.de/tee-care



	T-W07401	T-T08436	T-T08466	T-M07432	T-T07432	T-T07456	T-T07466	T-T07476
Enzyme pre-cleaning			✓				✓	✓
Immersion disinfection		✓	✓		✓	✓	✓	✓
Leak test				✓		✓		✓
Storage	✓	✓	✓					



WITH MOUNTING OPTION FOR PLUG PROTECTION





see page 15

Accessories for mobile units





Z-Z00300



Z-Z00310



Z-Z00320



Z-Z00330





T1 Z-Z00020



T2 Z-Z00030



T4 Z-Z00040



TU Z-Z00050



TU with optional GE plug protection type D $\,$ Z-Z00060




 $\textbf{Special version on request} \quad \text{Z-Z00070}$





Z-Z00100



Z-Z00110



Z-Z00120



Z-Z00130



see page 5







Z-TW0005

SCHOOL HOLDER



Z-TW0002

DISINFECTANTDISPENSER



Z-TW0001









HOLDERS

for test equipment

Z-Z00210





80 mm × 75 mm × 115 mm (L × w × H) Z-Z00230





Z-TT0006

DISCARD/SECRETION BOWL



Z-TW0003

NOTICE DRAIN VALVE



Z-TT0005

NOTICE PROTECTION SYSTEMS



GE Z-Z00007



Philips Z-Z00008



Special version on request Z-Z00009

Accessories for stationary units









Z-Z00630



7-700640



Z-Z00650



Z-Z00670







for GE Bite Hole Indicator Z-Z00220



Z-Z00240

N PROTECTORS

for stationary "DIRECT+" hygiene rear panels



T-W73416




T-W73426





Plug holders are available for the following manufacturers:



- + Philips
- + GE
- + Siemens
- + Canon
- + Fujifilm Sonosite
- + Mindray
- + Aloka Hitachi
- + BK Medical
- + Esaote
- + Samsung
- + Toshiba

Please contact us if your probe is not listed – we can develop a custom plug holder for you.



17

Easy online configuration: https://configurator.provita.de/tee-care

OVERVIEW

Dimensions of the different plug holders


	Designation	Standard rail 25 mm × 10 mm	Protector & hygiene rear panel
85 85	Plug holder T1	Z-Z00020	Z-Z00620
85 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	Plug holder T2	Z-Z00030	Z-Z00630
100	Plug holder T4	Z-Z00040	Z-Z00640
112,5 135.1	Plug holder GE/TU (1)	Z-Z00050	Z-Z00650
112,5	Plug holder GE/TU (2)	Z-Z00060	

18









Do you know our other catalogues?



provita medical gmbh & co. kg
Auf der Huhfuhr 8
42929 Wermelskirchen
Germany

Phone +49 (0) 21 93 / 51 05 - 0 Fax +49 (0) 21 93 / 51 05 - 269 info@provita.de www.provita.de



