

# ENDO Protector for standard rail

**500 mm**

**E-C11109**

## PRODUCT DESCRIPTION

<b>Handle bearing (s)</b>	2 units, 1 x Z-HA0001, 1 x Z-HA0005
<b>Working length</b>	500 mm
<b>Disinfection</b>	Thermal disinfection up to 95°C
<b>Material of housing</b>	Aluminum with special coating for thermidisinfection up to 95°C
<b>Endo Protect</b>	Storage options for endoscopes
<b>Marking insert</b>	white = neutral
<b>Type of Fixation</b>	Mobile, 2 standard rail adapters for equipment rails 25 x 10 mm
<b>Total load capacity (kg)</b>	2

## PACKAGING UNIT: UNIT

Width Box1	0,20 m
Height Box1	0,20 m
Length Box1	0,60 m
Weight Box1	2,00 kg



## ENDO Protector for standard rail



### Area of application

Dust protected and mechanically protected storage and transport of all endoscopes

Immediate usability of the endoscope

Mobile and stationary use



### Technology

Mechanically protected storage

Customizable through accessories



### Safety

Optics / Distal end of scope protected against impact and damage

Smooth surfaces for a safe disinfection result

Thermo-disinfectable  
Protector enables a validatable disinfection result

Thermo-disinfection at 95 °C of the protector tested with 1000 reprocessing cycles