

# IV-Pole Intensive care

## Material: stainless steel

I-I13212

### PRODUCT DESCRIPTION

|                               |  |
|-------------------------------|--|
| Field of application          | Intensive Care                             |
| Material                      | Stainless steel                            |
| Colour of plastic parts       | Blue                                       |
| Bottle Holder                 | Stainless steel to screw in, bottle holder |
| Number of hooks               | 4  |
| Arrangement of the hooks      | 90° to each other                          |
| Type of height adjustment     | One-hand height adjustment                 |
| Maximum length/height         | 2,350                                      |
| Minimum length/height         | 1,500                                      |
| Inner tube diameter           | Ø 18 mm                                    |
| Outer tube diameter           | Ø 25 mm                                    |
| Length of outer tube          | 1,200                                      |
| Base leg profile              | 45 x 25 mm                                 |
| Base weight                   | Weighted                                   |
| Base diameter                 | Ø 680 mm                                   |
| Castors                       | Apparatus castors                          |
| Castor diameter               | Ø 75 mm                                    |
| Number of castors             | 5; 2 with brakes and elect. conductivity   |
| Drip vessel sterilizable      | Yes  |
| Total load capacity (kg)      | 21   |
| Load capacity per system (kg) | 8  |
| Load capacity per hook (kg)   | 2  |



### PACKAGING UNIT: UNIT

|             |        |
|-------------|--------|
| Width Box1  | 0,22 m |
| Height Box1 | 0,12 m |
| Length Box1 | 1,47 m |
| Weight Box1 | 9,1 kg |



Best hygienic properties  
stainless steel (1.4301)



I-I13212

## IV-Pole Intensive care



### Economy

Compact packaging allows for low-cost shipping



### Safety

The base legs are sealed with plastic plugs (no risk of injury)

Low risk of pinching when retracting and extending the IV-pole thanks to the built-in safety distance

an integrated base-weights facilitates a safer attachment of various pumps and other medical equipment

No risk of injury thanks to the rounded shape of the bottle hooks

The geometry of the hooks prevents the infusion bottle from falling out

