

# IV Stand Normal care, heavy weighted Material: chromed steel

# I-N41222

### **PRODUCT DESCRIPTION**

Normal Care
Chromed steel
Blue
Chrome-plated steel bottle cross
4
90° to each other
Screw height adjustment
2.090 mm ± 10mm
1.250 mm ± 10mm
Ø 18 mm
Ø 25 mm
1.000 mm
25 x 25 mm
Weighted
Ø 635 mm
Twin castors
Ø 50 mm
5; 2 with brakes and elect. conductivity
Yes
11
8
2



#### **PACKAGING UNIT: UNIT**

Width Box1	0,19 m
Height Box1	0,10 m
Length Box1	1,19 m
Weight Box1	6,80 kg





many different configurations



© provita medical gmbh & co. kg – Technische Änderungen vorbehalten. https://www.provita.de/detail/5f872bf8a4f94cd58320d08b6212eb51

# I-N41222



# IV Stand Normal care, heavy weighted









optimal hygienic characteristics (wipe disinfection)

very easy and fast assembly through a plug in and out system with only 1 screw connection

comfortable height adjustment through a screw-height adjustment with an integrated, surface-gently mold clamping element



All plastic parts are lightfast

castors with non-marking tread

The integrated pull-out limiter prevents the adjustable tubes from separating

Safety bottle cross according to DIN ISO 15375

a sealing cap prevents the IVpole from hygienical senstive liquids



Compact packaging allows for low-cost shipping



The base legs are sealed with plastic plugs (no risk of injury)

Low risk of pinching when retracting and extending the IV-pole thanks to the built-in safety distance

an integrated base-weights facilitates a safer attachment of various pumps and other medical equipment

No risk of injury thanks to the rounded shape of the bottle hooks

The geometry of the hooks prevents the infusion bottle from falling out

