

# Mobile unit "MOBILE+" with 2 tubes

## Pre-cleaning and disinfection by immersion

T-T05060

### PRODUCT DESCRIPTION


Cleaning/disinfection tubes	Cleaning and Immersion disinfection
Field of application	TEE-Care (Ultrasound probes)
Dip tube length	1.160 (ca.)
Tube diameter	Ø 33,7 mm
Material support tube	stainless steel
Length of supp. tube	1.650
Function of the unit	disinfection including pre-cleaning
Reprocessing of the probe	by immersion disinfection
Base diameter	Ø 630 mm
Castors	Twin castors
Total load capacity (kg)	18

### PACKAGING UNIT: UNIT

Width Box1	0,7 m
Height Box1	0,2 m
Length Box1	0,7 m
Weight Box1	10 kg
Width Box2	0,3 m
Height Box2	0,3 m
Length Box2	1,7 m
Weight Box2	6 kg



Time-saving cleaning



T-T05060

## Mobile unit "MOBILE+" with 2 tubes

